




Macherey & Nagel QUANTOFIX sulfite

Product Images



Description

QUANTOFIX Sulfite test strips 6 x 95 mm Measuring range: 0-10-25-50-100-250- 500-1000 mg/L SO₃²⁻ Pg. à 100 determinations

Additional Information

No.	ML-6319
Dual use indicator	false
Manufacturer (Brand)	Macherey & Nagel
VGKL number	244691306
EAN	4058072965075
Gross weight	0.042
Manufacturer EAN	4046681004643
Storage temperature	Raumtemperatur
Net weight	0.04200 kg
S-phrases	S 22
Country of origin	Deutschland
Region of origin	Nordrhein-Westfalen
Packing width	0.04000 m
sterile	No
Packing height	0.04000 m
Packing depth	0.12300 m
Packaging volume	0.00019700 ccm
Customs tariff number	38229000

