



Macherey & Nagel QUANTOFIX Zinc

Product Images



Description

QUANTOFIX Zinc test strips 6 x 95 mm Measuring range: 0-2-5-10-25-50-100 mg/L Zn²⁺ Pg. à 100 determinations incl. reagents and reaction vessel __UN 3316 Chemistry test kit 9 II 0,060 L ADR/GGVSEB M11 (E) ADR 3.3.1/251: LQ PAX 960 10 Kg CAO 960 10 Kg

Additional Information

| | |
|-----------------------|--|
| No. | ML-6322 |
| Dual use indicator | false |
| Manufacturer (Brand) | Macherey & Nagel |
| VGKL number | 244691310 |
| GHS symbol | GHS05 |
| EAN | 4058072965105 |
| Gross weight | 0.228 |
| Manufacturer EAN | 4046681004681 |
| H-phrases | H314 |
| Storage temperature | Raumtemperatur |
| Net weight | 0.22800 kg |
| P phrases | P260sh. P280sh. P303. 361. 353. P305. 351. 338. P310 |
| R-phrases | R 35 |
| S-phrases | S 22-26-37/39-45 |
| UN number | 3316 |
| Country of origin | Deutschland |
| Region of origin | Nordrhein-Westfalen |
| Packing width | 0.17000 m |
| sterile | No |
| Packing height | 0.04800 m |
| Packing depth | 0.12500 m |
| Packaging volume | 0.00102000 ccm |
| Customs tariff number | 38229000 |
| Signal word | Gefahr |

