



Macherey & Nail rDNase Set

Product Images

Description

rDNase Set Recombinant DNase and rDNase Reaction Buffer, For RNA minipreparations with NucleoSpin kits no dangerous goods according to transport regulations

Additional Information

No.	ML-4498
Dual use indicator	false
Manufacturer (Brand)	Macherey & Nagel
EAN	4058072946760
Gross weight	0.1046
H-phrases	H334
Storage temperature	Raumtemperatur
Net weight	0.10460 kg
P phrases	P261sh. P342. 311
R-phrases	R 42/43
S-phrases	S 22
Country of origin	Deutschland
Region of origin	Nordrhein-Westfalen
Packing width	0.15000 m
sterile	No
Packing height	0.08500 m
Packing depth	0.11000 m
Packaging volume	0.00140300 ccm
Customs tariff number	38229000

