



Manitowoc Storage D-320

Product Images



Description

Storage unit suitable for modular flake ice maker RFP 0320A & RFP 0620A; for optimal utilization of required ice capacities Housing: DuraTech coating/ Capacity: 95kg Height adjustable feet Weight: 34kg Dimensions W x D x H: 559 x 864 x 965mm

Additional Information

No.	TW-0003
Manufacturer (Brand)	Manitowoc
VGKL number	101702320
EAN	4058072972844
Color	silber
Material	Kunststoff (unspezifiziert)
Wide	559 mm
Height	965 mm
Depth	864 mm

