



Manitowoc Storage D-420

Product Images



Description

Storage unit suitable for modular flake ice maker RFP 0320A & RFP 0620A; for optimal utilization of required ice capacities Housing: DuraTech coating/ Capacity: 141kg Height adjustable feet Weight: 49kg Dimensions W x D x H: 559 x 864 x 1270mm

Additional Information

No.	TW-0004
Manufacturer (Brand)	Manitowoc
VGKL number	101702420
EAN	4058072972851
Color	silber
Material	Kunststoff (unspezifiziert)
Wide	559 mm
Height	1270 mm
Depth	864 mm

