



Manitowoc Storage D-570

Product Images



Description

Storage unit suitable for the modular flake ice maker RFP 1300 for optimal utilization of the required ice capacities Housing: DuraTech coating/ Capacity: 195 kg Height-adjustable feet Weight: 52 kg Dimensions W x D x H: 762 x 864 x 1270mm

Additional Information

No.	TW-0005
Manufacturer (Brand)	Manitowoc
EAN	4058072972868
Color	silber
Material	Kunststoff (unspezifiziert)
Wide	762 mm
Height	1270 mm
Depth	864 mm

