



Manitowoc Storage D-970

Product Images



Description

Accumulator suitable for modular flake ice maker RFP 1300* * with adapter K-00470 for optimal utilization of required ice capacities
Housing: DuraTech coating/ Capacity: 322 kg Height adjustable feet Weight: 68 kg Dimensions W x D x H: 1,219 x 864 x 1270mm

Additional Information

No.	TW-0006
Manufacturer (Brand)	Manitowoc
EAN	4058072972875
Color	silber
Material	Kunststoff (unspezifiziert)
Wide	1219 mm
Height	1270 mm
Depth	864 mm

