



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex screw regulator for volume adjustment for Calibrex™ organo 525 dispenser, yellow

€14.04
plus VAT &
Shipping

Product Images



Description

Screw knob for adjusting the volume

As an alternative adjustment knob to the preassembled slide regulator, depending on user preference. The preassembled slide control can be exchanged for a screw control.

Additional Information

| | |
|------------------------|--------------------|
| No. | SR-0251 |
| Manufacturer (Brand) | Socorex |
| old neolab article no. | 7-7630 |
| EAN | 4058072367695 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | No |
| Color | yellow |
| Material | Polypropylene (PP) |
| sterile | No |
| suitable for | Socorex Dispenser |
| Manufacturer | Socorex |

