



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex set for increased piston seal (piston + cylinder), for 5 ml Dosys™ self-filling syringes

€47.63
plus VAT &
Shipping

Product Images



Description

Set for increased piston tightness

A special glass cylinder and piston with X-ring increase the syringe tightness if the application requires it (e.g. filtration, aspiration, etc.)

(e.g. filtration, aspiration, etc.).

Additional Information

No.	SR-0300
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7701
EAN	4058072368173
Transport temperature	Room temperature
sterile	No
suitable for	Socorex self-filling syringes
Manufacturer	Socorex

