



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex replacement capillary tubes for Socorex Acura® 846 direct displacement micropipettes 1 - 5 µL, 1xwhite

€69.99
plus VAT &
Shipping

Product Images



Description

Replacement glass capillaries for Socorex direct displacement micropipette Acura®capillar 846

The capillaries have color coding that matches the corresponding micropipette.
Supplied in a convenient box with lid.

Additional Information

No.	SR-0137
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7404
EAN	4058072366674
Transport temperature	Room temperature
sterile	No
Length	100 mm
suitable for	Socorex pipettes
Manufacturer	Socorex

