



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex replacement capillary tubes for Socorex Acura® 846 direct displacement micropipettes 10 – 50 µL, green

€69.99
plus VAT &
Shipping

Product Images



Description

Replacement glass capillaries for Socorex direct displacement micropipette Acura®capillar 846

The capillaries have color coding that matches the corresponding micropipette.
Supplied in a convenient box with lid.

Additional Information

| | |
|------------------------|------------------|
| No. | SR-0139 |
| Manufacturer (Brand) | Socorex |
| old neoLab article no. | 7-7406 |
| EAN | 4058072366698 |
| Transport temperature | Room temperature |
| sterile | No |
| Length | 100 mm |
| suitable for | Socorex pipettes |
| Manufacturer | Socorex |

