



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Socorex replacement capillary tubes for Socorex Acura® 846 direct displacement micropipettes 60 - 100 µL, blue

**€140.09**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### Replacement glass capillaries for Socorex direct displacement micropipette Acura®capillar 846

The capillaries have color coding that matches the corresponding micropipette.  
Supplied in a convenient box with lid.

## Additional Information

---

No.	SR-0140
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7407
EAN	4058072366704
Transport temperature	Room temperature
sterile	No
Length	100 mm
suitable for	Socorex pipettes
Manufacturer	Socorex

