



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

**Socorex replacement flask and glass capillary  
for Socorex Acura® 846 direct displacement  
micropipettes 5 - 25 µL, 2xwhite**

**€57.72  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

### Replacement flask with glass capillaries

Replacement flask with matching glass capillary for <sup>Acura®</sup> 846 direct displacement micropipettor. Each flask is supplied with a glass capillary.

The capillaries and flasks are color coded to match the corresponding micropipette Micropipette.

## Additional Information

---

No.	SR-0143
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7410
EAN	4058072366735
Transport temperature	Room temperature
sterile	No
suitable for	Socorex pipettes
Manufacturer	Socorex

