



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Socorex replacement flask and glass capillary
for Socorex Acura® 846 direct displacement
micropipettes 100 – 200 µL, red**

**€57.72
plus VAT &
Shipping**

Product Images



Description

Replacement flask with glass capillaries

Replacement flask with matching glass capillary for ^{Acura®} 846 direct displacement micropipettor. Each flask is supplied with a glass capillary.

The capillaries and flasks are color coded to match the corresponding micropipette Micropipette.

Additional Information

No.	SR-0146
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7413
EAN	4058072366766
Transport temperature	Room temperature
sterile	No
suitable for	Socorex pipettes
Manufacturer	Socorex

