



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex Suction Tube Conversion Set for Dosys™ Self-Filling Syringes

€29.95
plus VAT &
Shipping

Product Images



Description

Bottle holder conversion set for alternative liquid supply

Set includes bottle holder, hose attachment, feed needle, blocking ring, rubber seal (improves density between bottle and holder when bottle neck is thinner than standard) and slider.

No tools required.

Additional Information

No.	SR-0299
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7700
EAN	4058072368166
Transport temperature	Room temperature
sterile	No
Manufacturer	Socorex

