



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex volume unit for Acura® electro reduced shaft length, 50µl, blue only 927

€63.86
plus VAT &
Shipping

Product Images



Description

Socorex volume units for single-channel microliter pipettes Acura® electro 926XS Reduced shaft length

The control element of Socorex electronic pipettes fits a choice of 27 different volume units, available separately. This extends the working range and makes electronic pipetting accessible to all budgets.

Calibration specifications of additional units – according to factory control certificate or by the pipette owner – remain stored in the microprocessor.

Additional Information

No.	SR-0060
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7256
EAN	4058072366117
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	light gray
Material	Polypropylene (PP)
Volume MAX	50 µl
Volume MIN	2.5 µl
sterile	No
suitable for	Socorex pipettes
Number of channels	1
Manufacturer	Socorex
Model Microliter pipette/ Dispenser	Socorex Acura® electro XS 926

