



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Socorex volume unit for Acura® electro standard length, 10µl yellow, blue only 927

**€63.86**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### Socorex volume units for single-channel microliter pipettes Acura® electro 926XS Regular shaft length

The control element of Socorex electronic pipettes fits a choice of 27 different volume units, available separately. This expands the working range and makes electronic pipetting accessible accessible to all budgets.

Calibration specifications of additional units – according to factory control certificate or by the pipette owner – are stored in the microprocessor.

## Additional Information

---

No.	SR-0057
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7253
EAN	4058072366087
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	light gray
Material	Polypropylene (PP)
Volume MAX	10 µl
Volume MIN	0.5 µL
sterile	No
suitable for	Socorex pipettes
Number of channels	1
Manufacturer	Socorex
Model Microliter pipette/ Dispenser	Socorex Acura® electro XS 926

