




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex volume unit for Acura® electro 10 ml, blue 937

€86.42
plus VAT &
Shipping

Product Images



Description

Socorex volume units for Acura® electro 936 single-channel macroliter pipettes

The control element of the Socorex electronic pipettes fits a choice of 27 different volume units, available separately. This expands the working range and makes electronic pipetting accessible accessible to all budgets.

Calibration specifications of additional units - according to factory control certificate or by the pipette owner - are stored in the microprocessor.

Additional Information

No.	SR-0070
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7266
EAN	4058072366216
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	light gray
Material	Polypropylene (PP)
Volume MAX	10 ml
Volume MIN	0.5 ml
sterile	No
suitable for	Socorex pipettes
Number of channels	1
Model Microliter pipette/ Dispenser	Socorex Acura® electro XS 936
Manufacturer	Socorex

