



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® replacement rubbers for 15 ml tubes

€2.51
plus VAT &
Shipping

Product Images



Description

Replacement rubber holder for the racks of the Sunlab Diggi Mixer (D-8801).

Suitable for 15 ml tubes.

By using the rubber rings, the reaction tubes and centrifuge tubes in different sizes can be better fixed in the rack.

Additional Information

No.	D-8820
Manufacturer (Brand)	neoLab
EAN	4058072159481
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	Black
Material	Rubber
sterile	No
suitable for	Test tubes/tubes
Manufacturer	Sunlab

