



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® replacement rubbers for 50 ml tubes

**€2.51**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Replacement rubber holder for the racks of the Sunlab Diggi Mixer (D-8801).

Suitable for 50 ml tubes.

By using the rubber rings, the reaction tubes and centrifuge tubes in different sizes can be better fixed in the rack.

## Additional Information

---

No.	D-8821
Manufacturer (Brand)	neoLab
EAN	4058072158132
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	Black
Material	Rubber
sterile	No
suitable for	Test tubes/tubes
Manufacturer	Sunlab

