



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657




neoLab® replacement rubbers for 1.5 ml tubes

**Price on
request**

**plus VAT &
Shipping**

Product Images



Description

Replacement rubber holder for the racks of the Sunlab Diggi Mixer (D-8801).

Fits 1.5 ml.

By using the rubber rings, the reaction tubes and centrifuge tubes in different sizes can be better fixed in the rack.

Additional Information

No.	D-8822
Manufacturer (Brand)	neoLab
EAN	4058072156930
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	Black
Material	Rubber
sterile	No
suitable for	Test tubes/tubes
Manufacturer	Sunlab

