



# Systec steam/air mixture process

## Product Images

## Description

for sterilization of liquids in closed, pressure-sensitive containers. Prevents deformation of the container to high by internal pressure during sterilization. Rapid re-cooling not included. For Systec D-Series from 65 liters\* and Systec V-Series from 65 liters\* \*below 65 liters on request; For model series: VX/DX

## Additional Information

|                      |               |
|----------------------|---------------|
| No.                  | SX-0042       |
| Manufacturer (Brand) | Systec        |
| EAN                  | 4066292306653 |

