



Stakpure Ion Exchanger Mini, pressureless

Product Images



Description

The ideal water deionizer for any laboratory where small quantities of water are required.

Pressureless ion exchange system filled with mixed bed resin.

The unit is supplied complete with an analog conductivity meter, wall bracket and hose set.

Suitable for small laboratories with a daily requirement of up to 10 liters.

Ion exchange systems comply with the **DIN EN 285** standard when used in the upstream part of the water circuit for use in sterilizers.

Filter system protects against resin leakage.

Mini 450 model incl. cartridge.

Additional Information

No.	PE-0070
Manufacturer (Brand)	Stakpure
old neoLab article no.	2-0147
EAN	4066292550421
diameter	100 mm
Height	600 mm
Weight	3 kg
Connection water	3/4"
Application product water	Allgemeine Laboranwendungen (z.B. Autoklavieren. Medienzubereitung. etc.)
Feed water pressure	1 bis 6 bar
Spare part	PE-0071
Filling cartridge	Mischbettharz (Kationen-/Anionenharz)
Capacity at 10° GSG	425 l
Marking by manufacturer with	ASTM
Conductivity deionized water	0.1 bis 20 µS/cm
Model Ion exchanger	neoLab Ionentauscher Mini (2-0147)
Production capacity deionized water	50 l/h
Quality product water	VE-Wasser (Qualität abhängig vom Speisewasser)

