



Thermo Scientific Cryoboxes

Product Images

Description

Cryoboxes Meet every need and budget with Thermo Scientific™ cryoboxes. Description: Polycarbonate, 132 x 132 x 52 mm
Construction/Design: Nunc MAX-100 Polycarbonate CryoStore Boxes Shape: Square Type: cryobox Capacity: 100 (1.0 to 1.8 ml vials)
Quantity: 4 pieces/pack; 24 pieces/carton Tube Array: 10 x 10 Dimensions (L x W x H): 132 x 132 x 52 mm (5.25 x 5.25 x 2.13 in.)
Length (anglo): 5.25 in. Length (metric): 132 mm Width (Anglo): 5.25 in. Width (metric): 132 mm Height (Anglo): 2.13 in. Height (metric): 52 mm Unit Size: Case of 24

Additional Information

No.	FE-0723
Manufacturer (Brand)	Thermo Scientific
VGKL number	2900374187
EAN	4066292514348
Gross weight	4.11
Net weight	1.00000 kg
Color	grau
Material	Polycarbonat (PC)
Country of origin	Mexiko
Packaging volume	0.02600000 ccm
Wide	132 mm
Customs tariff number	39269097
Height	52 mm
Length	132 mm

