



Thermo Scientific Reverse Osmosis Membrane

Product Images

Description

Reverse Osmosis Membrane Select the required Thermo Scientific™ Reverse Osmosis (RO) membrane for your Thermo Scientific™ Barnstead™ water treatment system. Description: RO membrane with integrated pretreatment module, 6 l/h For use with (devices): Smart2Pure Type: Reverse Osmosis Membrane Unit Size: Each

Additional Information

No.	LD-0065
Manufacturer (Brand)	Thermo Scientific
VGKL number	108662006
EAN	4066292000599

