



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine level sensor, for Maxi aspiration device

€227.40
plus VAT &
Shipping

Product Images



Description

Capacitive level sensor that switches off the pump when the preset liquid level is reached (the sensor is mounted outside the bottle and is therefore not contaminated).

including screwdriver 1-7230

Additional Information

No.	2-9332
Manufacturer (Brand)	neoLabLine
EAN	4058072070052
Transport temperature	Room temperature

