



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



IST Innuscreen innuPREP DNase I Digest Kit, 250 reactions

€388.50
plus VAT &
Shipping

Product Images



Description

The innuPREP DNase I Digest Kit is used for the highly efficient removal of DNA from appropriately contaminated RNA samples. Enzymatic DNA digestion is designed to be conveniently performed during RNA isolation.

- Can be quickly and easily integrated into the RNA isolation process
- Easy handling due to direct DNase I digestion on the column
- Excellent purity of eluted RNA due to efficient DNA digestion
- Free of RNase and other interfering factors

Additional Information

No.	C-6155
Manufacturer (Brand)	IST Innuscreen
EAN	4058072054960
Transport temperature	Wet ice

