



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neoFroxx XXprep Kit for DNA fragments (from PCR reactions & agarose gels) 250 reactions

€432.60 As low
as
€383.15
plus VAT &
Shipping

Product Images



Description

The **XXprep Kit for DNA fragments** (from PCR reactions & agarose gels) has been developed for the efficient extraction of DNA fragments from TAE or TBE agarose gels and utilizes a novel 2-step technology for purification of amplification products from PCR reaction mixtures and sequencing reaction mixtures.

The Spin Filter column has a binding capacity of over 20 µg, allowing high yield and excellent quality in DNA fragment isolation.

Recovery rates are up to 95% depending on fragment size.

Other features:

easy handling, flexible choice of elution volume and fast purification protocols

Advantages:

- Combi Kit
- PCR purification in a two-step procedure (without unnecessary washing)
- reduced use of plastics
- minimal elution volumes (up to 10 µl)

Additional Information

No.	5482RC250
Manufacturer (Brand)	neoFroxx
EAN	4066306002403
Transport temperature	Room temperature
autoclavable at 121°C	No
sterile	No

