



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® EZ Dry desiccant bag for 100 ml, 5 pcs./pack

Special Price

€35.00 was

€52.50

**plus VAT &
Shipping**

Product Images



Description

The neoLab-EZ Dry desiccant bags filled with molecular sieve are ideally suited for the preparation and maintenance of water-free solvents. Anhydrous conditions are ensured from day one where it is needed (e.g. oligonucleotide synthesis). After use, the bag can be easily removed with tweezers. The sometimes complex filtration with molecular sieve is thus bypassed

- **For volume: 100 ml**
- **For smallest bottle neck: 12 mm**

Additional Information

No.	1-6420
Manufacturer (Brand)	neoLab
EAN	4058072087272
autoclavable at 121°C	No

