



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Drying beads w. color indicator, 1 l

€30.00
plus VAT &
Shipping

Product Images



Description

Indicator silica gel adsorbent ^{KC-Trockenperlen®} for moisture control.

Spherical silica gel with very high mechanical strength. Free of heavy metals and environmentally compatible.

A very well visible color change guarantees a clear moisture indication. The degree of moisture absorption is indicated by a color change from orange to colorless at approx. 6 wt.% water content. Even after this, the drying beads continue to adsorb water until the maximum absorption capacity is reached.

They can be dried again in the drying oven at 100° to 110°C and can then be used again.

Grain size: 2 to 5 mm

Additional Information

No.	1-7031
Manufacturer (Brand)	neoLabLine
old neoLab article no.	1-7031
EAN	4066292697287

