



CellCamper® Cooling Unit Cool (0°-4°C)

Product Images



Description

The cooling units for the CellCamper® Midi and Maxi cover three different temperature ranges (0 to +4°C / -20 to 0°C / -22 to -18°C), which can be maintained between 6 and 14 hours.

Cooling unit **Cool (2-3716)** :

- Temperature range: 0 to +4°C
- Duration: up to 14 hours
- Cooling time at -20°C: 4 hours
- Cooling time at -80°C: 2 hours

After cooling at -80°C, leave at room temperature for 10 minutes to eliminate surface temperature.

Additional Information

No.	2-3716
Manufacturer (Brand)	CellCamper
EAN	4058072090043

