




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



CellCamper® Cooling Unit Frozen (-20°-0°C)


Special Price

€179.00 was

€239.00

**plus VAT &
Shipping**

Product Images



Description

The cooling units for the CellCamper® Midi and Maxi cover three different temperature ranges (0 to +4°C / -20 to 0°C / -22 to -18°C), which can be maintained between 6 and 14 hours.

Cooling unit **Frozen (2-3717)** :

- Temperature range: -20 to 0°C
- Duration: up to 12 hours
- Cooling time at -20°C: 7 hours
- Cooling time at -80°C: 3.5 hours

After cooling at -80°C, leave at room temperature for 10 minutes to eliminate surface temperature.

Additional Information

No.	2-3717
Manufacturer (Brand)	CellCamper
EAN	4058072090050
Transport temperature	Room temperature

