



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



CellCamper® Cooling Unit Ultrafrozen (-22°- -18°C)

Special Price

€194.00 was

€255.90

**plus VAT &
Shipping**

Product Images



Description

The refrigeration units for the CellCamper® Midi and Maxi cover three different temperature ranges (0 to +4°C / -20 to 0°C / -22 to -18°C) that can be maintained between 6 and 14 hours.

Ultrafrozen refrigeration unit(2-3718) :

- Temperature range: -22 to -18°C
 - Duration: up to 6 hours
 - Cooling time at -80°C: 4 hours
- After cooling at -80°C, leave at room temperature for 10 minutes to eliminate surface temperature.

Additional Information

No.	2-3718
Manufacturer (Brand)	CellCamper
EAN	4058072090067
Transport temperature	Room temperature

