



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN double gooseneck lighting (UK version) 6W LED: 5600-6300 K: 300mm

€165.00
plus VAT &
Shipping

Product Images



Description

- To offer you maximum flexibility and comfort in stereomicroscopy, select your favorite illumination here
- These professional illumination units guarantee outstanding light quality with constant intensity on the object
- Whether space-saving ring illumination or cold light sources with light guides, our range leaves nothing to be desired
- With the polarization ring illumination OZB-A7101, you also have an outstanding component at your disposal, especially optimized for viewing glossy surfaces
- Of course, these external illumination units are also suitable for your standard stereo microscope (exception: the ring illuminators do not fit in combination with the following series: OSE-1, OSF-4G, OZL-44, OZL-45, OZL-45R, OZC-5 and OZG-4)
- NEW: The OZB-A4516 LED gooseneck illuminator with focusable light cone allows you to customize your illumination. A punctual or a scattered beam allows you the optimal illumination of your sample

Additional Information

No.	KP-5165
Manufacturer (Brand)	KERN
EAN	4058072147112
Transport temperature	Room temperature

