



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN microscope camera ODC,

€1,150.00
plus VAT &
Shipping

Product Images



Description

- The ODC-851 HDMI microscope camera is specially designed for direct HDMI connection to your HDMI-enabled playback device. Images can also be saved directly to the included SD card.
- As an alternative, a USB 2.0 port is also available for connecting a laptop or PC via the USB interface.
- Power is supplied via an external 12 V power unit.
- **Scope of delivery:** Camera, software, USB mouse, USB 2.0 cable (length: 2 m), HDMI cable (length: 2 m) and SD card (16 GB)
- Please order the C-mount adapter suitable for your KERN microscope at the same time.

Additional Information

No.	KP-5025
Manufacturer (Brand)	KERN
VGKL number	105830851
EAN	4058072145804
Transport temperature	Room temperature

