



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN stereo zoom microscope head OZO, 0.800 x – 7 x

**Price on  
request  
plus VAT &  
Shipping**

### Product Images



## Description

- To give you full flexibility for your special needs and applications, we offer here a wide range of stereo microscope heads, universal stands and external illuminators, which are easy to combine.
- Due to the different features of the stereo microscope heads, as well as the flexibility of the universal stands and the professional fixation of our mounts, you can configure your desired microscope as you like.
- For this purpose, various microscope heads from our product lines are available to you, each as a binocular or trinocular version.
- For the connection of a camera to the trinocular version, a C-mount adapter is required, which has to be selected from the following model equipment list

**Technical data:**

- Optical system: Greenough optics

## Additional Information

No.	KP-5301
Manufacturer (Brand)	KERN
EAN	4058072148324
Transport temperature	Room temperature

