



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN Stereo Zoom Microscope Head OZM, 0.7 x - 4 5

€1,280.00
plus VAT &
Shipping

Product Images



Description

- To give you full flexibility for your special needs and applications, we offer here a wide range of stereo microscope heads, universal stands and external illuminators, which are easy to combine.
- Due to the different features of the stereo microscope heads, as well as the flexibility of the universal stands and the professional fixation of our mounts, you can configure your desired microscope as you like.
- For this purpose, various microscope heads from our product lines are available to you, each as a binocular or trinocular version.
- For the connection of a camera to the trinocular version, a C-mount adapter is required, which has to be selected from the following model equipment list

Technical data:

- Optical system: Greenough optics

Additional Information

No.	KP-5273
Manufacturer (Brand)	KERN
EAN	4058072148072
Transport temperature	Room temperature

