




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Zeiss tripod, K MAT.

Special Price

**€531.90** was

**€591.00**

**plus VAT &  
Shipping**

## Product Images



## Description

---

**Stand K MAT** - stand base T310xW200xH35 mm - work surface T195xW160 mm - interface d=84 mm for tables - column 250 mm with drive and handle, lifting range 145 mm - Stemi support d=76 mm, load capacity 5 Kg, friction adjustable - 2x socket for IVI and lamp K LED - AL control element for on/off/dimming - ESD protection function: antistatic surface, grounding socket - built-in power supply 12V DC 24W/100...240V AC/50...60Hz Additionally included: - B/W plastic plate D=84x5 mm - Country-specific power cable EURO C8 - Dust cover

### Compatible microscopes

This item is compatible with the microscope(s) listed below. Please note: The list also includes microscopes which are no longer available in the current program.

Stemi 305  
Stemi 508  
Stemi 305 Demo  
Stemi 508 Demo  
Stemi 305 trino  
Stemi 508 trino  
Stemi 508 doc  
Stemi 305 cam

## Additional Information

---

|                       |                  |
|-----------------------|------------------|
| No.                   | CZ-0577          |
| Manufacturer (Brand)  | Zeiss            |
| EAN                   | 4066292500389    |
| Transport temperature | Room temperature |

