



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



2mag control unit for stirring systems 7-3089 to 7-3093

€1,385.00
plus VAT &
Shipping

Product Images



Description

Control unit for magnetic stirring systems for cell cultures

Powerful control unit for connection by cable to the stirring systems.

10-stage power setting.

Digital display for speed and power setting, SoftStart, QuickSet, saving of last function status when switched off.

Durable stainless steel housing, easy to clean.

Technical data

Stirring power: 4-40 watts (10 steps)

Weight: approx. 1.4 kg

Operating voltage: 100-240 V/50-60 Hz/1.5 A

Additional Information

No.	ZU-0085
Manufacturer (Brand)	2mag
old neoLab article no.	7-3093
EAN	4058072444518
Transport temperature	Room temperature

