



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Mini water deionizer including cartridge

€479.00
plus VAT &
Shipping

Product Images

Description

Mixed bed water deionizer Mini, unpressurized for the production of deionized water from drinking water

- Pure water quality down to $< 0.1 \mu\text{s}/\text{cm}$
- Desalination with quality mixed bed resins
- Capacity specification applies to: Feed water with approx. $300 \mu\text{s}/\text{cm}$ ($3.57 \text{ mmol}/\text{z}$, approx. 10°dh) and a loading of up to $20 \mu\text{s}/\text{cm}$ in the pure water
- The capacity decreases at: higher salt content in the feed water and/or at lower permissible conductance in the pure water
- All hydraulic connections: $r \frac{3}{4}$ "
- compact and simple system
- Easy installation
- Quick and easy cartridge change
- No waste water

Additional Information

No.	2-1171
Manufacturer (Brand)	neoLab
EAN	4058072189297
Transport temperature	Room temperature

