



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN weight set stainless steel fine turned, 1 mg - 10000000 mg, M1

**€1,350.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

- Test weight material: milligram weights aluminum/new silver, individual weights brass finely turned or stainless steel finely turned
- Milligram weight sets 1 mg-500 mg are supplied in a separate plastic case

## Additional Information

---

No.	KP-1457
Manufacturer (Brand)	KERN
EAN	4058072127831
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver
Material	Stainless steel
sterile	No
Weight max	10000000 mg
Weight min	1 mg
Manufacturer	CORE

