



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN weight set stainless steel fine turned, 1 mg - 5000000 mg, M1

€920.00
plus VAT &
Shipping

Product Images



Description

- Test weight material: milligram weights aluminum/new silver, individual weights brass finely turned or stainless steel finely turned
- Milligram weight sets 1 mg-500 mg are supplied in a separate plastic case

Additional Information

No.	KP-1969
Manufacturer (Brand)	KERN
EAN	4058072129323
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver
Material	Stainless steel
sterile	No
Weight max	5000000 mg
Weight min	1 mg
Manufacturer	CORE

