



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN set of weights brass fine turned, 1 g - 500 g, M1

**Price on  
request  
plus VAT &  
Shipping**

### Product Images



## Description

---

- Test weight material: milligram weights aluminum/new silver, individual weights brass finely turned or stainless steel finely turned
- Milligram weight sets 1 mg-500 mg are supplied in a separate plastic case

## Additional Information

---

No.	KP-2663
Manufacturer (Brand)	KERN
EAN	4058072131432
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	Gold
Material	Brass
sterile	No
Weight max	500 g
Weight min	1 g
Manufacturer	CORE

