



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN weight set in wooden block cast iron lacquered, 1 g - 5000 g, M3

€670.00
plus VAT &
Shipping

Product Images



Description

- Material test weight: brass finely turned and cast iron, lacquered
- Material case: wooden block

Additional Information

No.	KP-0461
Manufacturer (Brand)	KERN
VGKL number	645836298
EAN	4058072125783
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	Black
Material	Cast iron
sterile	No
Weight max	5000 g
Weight min	1 g
Manufacturer	CORE

