



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN weight set in wooden case stainless steel polished, 1 g - 2000 g, E2

**€2,340.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

- Material test weight: single weights stainless steel polished, milligram weights
- Material case: wood padded. Milligram weights (1 mg - 500 mg) in removable plastic case
- The table shows the smallest and the largest weight in each case
- **Note:** Shipping weights without a protective case can most likely result in damage to the weights that exceeds the respective allowable error limit. In the worst case, your weights may become unusable. Therefore, we offer you the appropriate protective and storage cases that prevent damage and keep your weights safe during their service life. Should you nevertheless wish your weights to be shipped without such a protective case, please confirm this to KERN in writing.

## Additional Information

---

No.	KP-2720
Manufacturer (Brand)	KERN
EAN	4058072131746
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver
Material	Stainless steel
sterile	No
Weight max	2000 g
Weight min	1 g
Manufacturer	CORE

