



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN weight set in wooden case brass fine turned, 1 mg - 100000 mg, M1

**Price on  
request  
plus VAT &  
Shipping**

### Product Images



## Description

---

- Material test weight: single weights brass finely turned or brass finely turned, milligram weights aluminum (1 mg-5 mg) / nickel silver (10 mg-500 mg)
- Material case: wood

## Additional Information

---

No.	KP-3270
Manufacturer (Brand)	KERN
EAN	4058072133559
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	Gold
Material	Brass
sterile	No
Weight max	100000 mg
Weight min	1 mg
Manufacturer	CORE

