



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN weight set in wooden case, polished stainless steel, F1 1 g - 50 g, F1

**€440.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

- Test weight material: milligram weights aluminum/new silver, individual weights polished stainless steel
- Milligram weight sets 1 mg-500 mg are supplied in a separate plastic case

## Additional Information

---

No.	KP-3377
Manufacturer (Brand)	KERN
EAN	4058072133856
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver
Material	Stainless steel
sterile	No
Weight max	50 g
Weight min	1 g
Manufacturer	CORE

