



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN weight set in wooden case stainless steel fine turned, 1 mg - 500000 mg, F2

€1,420.00
plus VAT &
Shipping

Product Images



Description

- Material test weight: single weights stainless steel fine turned, milligram weights aluminum (1 mg - 5 mg) / nickel silver (10 mg - 500 mg)
- Material case: wood

Additional Information

| | |
|-----------------------|------------------|
| No. | KP-3388 |
| Manufacturer (Brand) | KERN |
| EAN | 4058072133948 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| Color | silver |
| Material | Stainless steel |
| sterile | No |
| Weight max | 5000000 mg |
| Weight min | 1 mg |
| Manufacturer | CORE |

