



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN weight set in wooden case stainless steel fine turned, 1 g - 50 g, M1

€153.00
plus VAT &
Shipping

Product Images



Description

- Material test weight: single weights brass fine turned or brass fine turned, milligram weights aluminum (1 mg-5 mg) / nickel silver (10 mg-500 mg)
- Material case: wood

Additional Information

No.	KP-3544
Manufacturer (Brand)	KERN
EAN	4058072134723
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver
Material	Stainless steel
sterile	No
Weight max	50 g
Weight min	1 g
Manufacturer	CORE

