



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN weight set in wooden case, polished stainless steel, 1 mg - 1 kg, E1

**€3,200.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

- Milligram weight sets 1 mg–500 mg are supplied in a separate plastic case
- Material test weight: polished stainless steel

## Additional Information

---

|                       |                  |
|-----------------------|------------------|
| No.                   | KP-3808          |
| Manufacturer (Brand)  | KERN             |
| EAN                   | 4058072135430    |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes              |
| Color                 | silver           |
| Material              | Stainless steel  |
| sterile               | No               |
| Weight max            | 1000000 mg       |
| Weight min            | 1 mg             |
| Manufacturer          | CORE             |

